TROLLEY MAINTENANCE LOG

Trolley Serial #:						
_Intake HEPA Change:/ Exhaust HEPA Changed:/						
		Air Exhaust Pre-	Intake Pre-filter & Hoses		Trolly	HEPA Filters
Date	Initials	Filter ∆			Certified	
enter date/initials and $\sqrt{ m boxes}$						
		E2wks	E3mo	E12mo	E12mo	PRN

Trolley Air Exhaust Pre-Filter is changed every 2 weeks

With trolley running pre-filter is spritzed w/Sporicidin & bagged-out.
 Trolley exhaust area is wiped w/Sporicidin & new filter installed.
 Dirty exhaust filter bag is spritzed w/Sporicidin & removed for disposal

Trolley Rack Hoses are changed every 3 months

1.Hoses are removed & capped with glove – do not create aerosols. 2.Hoses are spritzed w/Sporicidin & removed to dirty CW for cleaning/decon

Trolley Air Intake Pre-Filter is changed/ sanitized every 12 months

Trolley Certification every 12 months

Trolley HEPA Filters Δ /Decon = 3 yrs/as needed - indicate intake (IT) or exhaust (EX) HEPA