

# Medium pressure Liquid Chromatography (MPLC)



## TELEDYNE-Isco Torrent combiflash (UV 200-760 nm):

- ❖ large scale from 33 to 300g with fraction collector
- ❖ Normal (Silica)

## TELEDYNE-Isco Rf MPLC (UV 200-360 nm):

- ❖ scale from 0.2 to 33 g with fraction collector
- ❖ Normal (Silica) , reverse phase (C18)

# High Performance Liquid Chromatography (HPLC)



## Shimadzu semi-preparative HPLC for sample fractionation and purification (UV, ELSD and RID):

- ❖ From 1 to 100 mg
- ❖ Normal (silica), reverse phase (C18, CN, ...)